



World Class Textile Laboratory Instruments

"Aadarsh Technologies"

is in the field of catering world class testing instruments to the diverse needs of various industries such as Textile, Paper, Paint, Pigments & plastic industry.

Our rich experience has enabled us to understand the customer requirements and the challenges faced by them. We are dedicated to providing our customers with technologically advanced testing machines and unmatched quality backed by technical support and after sales service; in short, world class products and services.

We are led by the team of industry technocrats whose vast industry experience and market knowledge has helped us in offering an accurate range of products to our client's at the most competitive prices. Our company also ensure to deliver the range at shortest delivery timings to our customers.

We offer total quality solutions at every stage of manufacturing.

Customer Satisfaction

We are committed to provide the equipment, knowledge & support to our clients in the way they require to enhance their levels of productivity & achieve high quality in their work. We provide every customer value for money.

The equipments we supply are used in various organizations including leading institutions and reputed exporters throughout India. We further looking ahead to increase our market reach and deal with the clients based in overseas market.

This will further provide a reputation to the name of our company in global markets.



We Offer

- Tailor Made Machinery
- Soft Flow Dyeing Machine
- Loop Steamer
- Continuous Dyeing Range
- Laboratory Winch
- Yarn heat setting machine



**CUT & FIT
GSM CUTTER**



**MARTINDALE ABRASION
TESTER**



GSM CUTTER



**DIGITAL CSP/TENSILE
STRENGTH TESTER**



**BUTTON SNAP PULL
TESTER**





IR DYEING MACHINE

Using infrared radiations for heating & can be customized as per requirement



ROTA DYER

Can be used to dye fabric or yarn up to 98 PC 05 Controller



PADDING MANGLE

Used in dyeing, printing and finishing processes can be supplied in horizontal/vertical models



COLOR MATCHING BOOTH

Supplied with 5 Illuminants with inner color neutral grey and all lamps are imported



DIGITAL AIR PERMEABILITY TESTER

Used to test the air permeability of various woven fabrics, including industries fabrics, non-woven and other breathable fabric.



CROCK METER

Used for checking dry and wet rubbing fastness




DIGITAL BURSTING STRENGTH TESTER

Ideal for testing wide range of paper and packaging material



FLAMMABILITY TESTER

Used for measuring Surface flammability of textile and plastics



DIGITAL WRAP REEL

Ideal for preparing hank or lea of accurate length to test the count and CSP



HOT AIR OVEN

Maximum Temperature 25 Available in different sizes 0 C.



PILLING TESTER

Used for checking pilling properties of textile material



ELMENDORF TEARING TESTER

Used for checking tearing strength of fabric by falling pendulum type



Address :

PAP R-339, Near Antony Garage
Rabale MIDC Area, Rabale,
Navi Mumbai - 400701

Regd. Office :

A-2301, Evergreen Heights,
Near Ozone Valley, Parsik Nagar,
Kalwa (W)Thane - 400605

Contact No :

+91 - 9223191421 / +91 - 8691007010

Email :

sales@aadarshtech.com /
info@aadarshtech.com

www.aadarshtech.com